

## Dead Time

In TCSPC the term dead time refers to the time the TCSPC system needs armed again after detecting an event. During the dead time the TCSPC system is blind.

If for example two photons are detected with the dead time of the TCSPC device the second photon will be lost. This leads to the so called [Pile-up Effect](#).

Note that not only TCSPC devices but also photon counting devices exhibit dead times.

Copyright of this document belongs to PicoQuant GmbH. No parts of it may be reproduced, translated or transferred to third parties without written permission of PicoQuant GmbH. All information given here is reliable to our best knowledge. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearances are subject to change without notice.



PicoQuant GmbH  
Rudower Chaussee 29 (IGZ)  
12489 Berlin  
Germany

P +49-(0)30-1208820-89  
F +49-(0)30-1208820-90  
info@picoquant.com  
www.picoquant.com